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# FLUXmedicare®

## Photodynamic Therapy System

User Manual for the FLUXMEDICARE® V2 System

**MDB TEXINOV®**  
56 route de Ferrossière  
38110 SAINT DIDIER DE LA TOUR  
France

Tél. : +33 (0)4 74 94 44 75  
Email : [medical@texinov.fr](mailto:medical@texinov.fr)

# Contents

1.	Introduction .....	33
	1.1. Description of the system.....	33
	1.2. Expected use.....	36
	1.3. Therapeutic indications .....	36
	1.4. Contraindications .....	36
	1.5. Undesirable effects.....	36
	1.6. Precautions for use.....	37
2.	Warnings, symbols .....	37
	2.1. Warnings and safety remarks.....	37
	2.2. Safety linked to the light device .....	38
	2.3. Product labelling.....	38
	2.4. Laser warning label.....	39
	2.5. Transport and storage markings .....	39
3.	Technical characteristics .....	40
	3.1. Description .....	40
	3.2. General information .....	40
	3.3. Electrical and optical performances .....	40
	3.4. Physical characteristics of the light source .....	41
	3.5. Ambient conditions .....	41
	3.6. Transport and storage conditions .....	41
	3.7. Inspection of the system .....	42
	3.8. Preparation of the source before treatment .....	42
	3.9. Connection of the optical fibres.....	43
4.	Protocol for use.....	43
	4.1. Preparation of the patient's skin.....	43
	4.2. Installation of the FLUXMEDICARE® V2 TEXLIGHT1-V2 system .....	43
	4.3. Putting the intermediate unit and harness in place.....	44
	4.4. Use of the light source.....	45
	4.4.1. Powering up the source.....	46
	4.4.2. Display screen.....	47
	4.4.3. Initialisation .....	47
	4.4.4. Starting treatment.....	47
	4.4.5. Pausing .....	48
	4.4.6. Stopping before the end of treatment .....	48
	4.5. End of treatment .....	49
	4.6. Installation of other FLUXMEDICARE® V2 accessories .....	49
5.	Error messages .....	50
6.	What to do in case of emergency .....	51

7.	Replacement of fuses.....	51
8.	Cleaning.....	52
	8.1. General.....	52
	8.2. Source.....	52
	8.3. Cap.....	53
	8.4. The stretch net.....	53
	8.5. Luminous bands.....	53
9.	Control of the source.....	53
10.	Recycling.....	54
11.	Electromagnetic compatibility.....	54
12.	Regulation.....	55
13.	Manufacturer's contact details.....	56

# Forward

Thank you for buying the FLUXMEDICARE® V2 photodynamic therapy system.

Please read this manual carefully before using the system and keep it for future reference. The content of this manual is based on information updated after compilation.

This manual of [FLUXMEDICARE® V2](https://www.fluxmedicare.fr/traitemement-pdt/) is available on the product website by following the link: <https://www.fluxmedicare.fr/traitemement-pdt/>

## 1. Introduction

### 1.1. Description of the system

The FLUXMEDICARE® V2 system is a light device intended to be used solely for photodynamic therapy (PDT) for the treatment of skin lesions.

The FLUXMEDICARE® V2 system is made up of a red light generator - called the light source - and a band of luminous textile transmitting the light.

The light source and the luminous band are connected.

The light source is connected to the mains electricity supply.

The light emitted has a peak wavelength of 638 nm.

The light dose received by the patient for the 2.5 hours' treatment is 12 J/cm<sup>2</sup>.

The system is made up of two kits sold separately:

- The **FLUXMEDICARE® V2 No.1 kit** - Ref. LIGHTSOURCE02 - is made up of the 3 following elements:
  - A light source
  - A harness (HARN2) for attachment to the torso to hold the intermediate unit. This can be sold separately.
  - A cap (CAP2) which can be sold separately.

Reference	Description	Photographs
<b>LIGHTSOURCE02</b>	1 x source 1 x harness 1 x cap	

The following elements are also sold separately:

Reference	Description	Photographs
<b>HARN2</b>	1 x harness	
<b>CAP2</b>	1 x cap	



- The **FLUXMEDICARE® V2 No.2 kit** - Ref. FLUXMED X-V2 - including the light device is made up of the following elements:
  - 12 stretch nets, 20 cm in length (NET05 and NET07) for single use only,
  - 1 batch of 10 pieces of “FLUXMEDICARE® Protect” protective cover (FMPOX) for single use,
  - one band of luminous textile covered with opaque blue fabric (sheath) with 1 cable to be connected to the source (TEXTLIGHTX-V2), associated with a specific count key to operated the source.

Reference	Description	Use	
<b>NET05</b>	6 nets: size No.5	Single use only	
<b>NET07</b>	6 nets: size No.7	Single use only	
<b>TEXTLIGHTX-V2</b> <b>X is a number from 1 to 6 corresponding to the model</b>	1 luminous textile (see the models in the following table) 1 x related count key	Reusable	

Reference	Description	Photographs
<b>TEXTLIGHT1-V2</b>	Bilateral sheath Illuminated surface: 18 cm x 20 cm 1 x count key	
<b>TEXTLIGHT2-V2</b>	Unilateral sheath Illuminated surface: 18 cm x 20 cm 1 x count key	

Reference	Description	Photographs
<b>TEXTLIGHT3-V2</b>	Small unilateral sheath Illuminated surface: 6 cm x 20 cm 1 x count key	
<b>TEXTLIGHT4-V2</b>	Small bilateral sheath Illuminated surface: 6 cm x 20 cm 1 x count key	
<b>TEXTLIGHT6-V2</b>	Bilateral temple sheath 2 illuminated surfaces: 6 cm x 9 cm* 1 x count key	

The following elements can be ordered separately:

Reference	Description	Use	Non-contractual photos
<b>NET05</b>	6 nets: size No.5	Single use only	
<b>NET07</b>	6 nets: size No.7	Single use only	
<b>FMP01</b>	“FLUXMEDICARE® Protect” protective cover in non-woven material for luminous textile TEXTLIGHT1-V2 - Single size	Single use only	
<b>FMP02</b>	“FLUXMEDICARE® Protect” protective cover in non-woven material for luminous textile TEXTLIGHT2-V2 - Single size	Single use only	
<b>FMP03</b>	“FLUXMEDICARE® Protect” protective cover in non-woven material for luminous textile TEXTLIGHT3-V2 - Single size	Single use only	
<b>FMP04</b>	“FLUXMEDICARE® Protect” protective cover in non-woven material for luminous textile TEXTLIGHT4-V2 - Single size	Single use only	
<b>FMP06</b>	“FLUXMEDICARE® Protect” protective cover in non-woven material for luminous textile TEXTLIGHT6-V2 - Single size	Single use only	

## 1.2. Expected use

---

The FLUXMEDICARE® V2 system is intended to be used by health professionals practising in doctors' surgeries, hospitals and clinics for PDT (photodynamic therapy) treatments.

The treatment should only be carried out by an authorised medical professional trained in photodynamic therapy.

For the PDT treatment, and before exposure to the light, the area of skin to be treated is covered by a cream containing 5-aminolevulinic acid or an ALA ester or methyl aminolevulinate (not supplied with the FLUXMEDICARE® V2 system) acting as a photosensitiser. These creams are only available on prescription.

## 1.3. Therapeutic indications / benefit

---

Photodynamic therapy is used to treat fine, non-pigmented precancerous skin lesions (actinic keratosis) on the face and scalp (Olsen Grade 1 to 2).

The target population is adults only.

This devices come with low level of pain and equivalent level of performance when compared to alternative treatment options and patient can move during the treatment.

## 1.4. Contraindications

---

The following contraindications are found for the phototherapy carried out in association with a pharmaceutical cream containing 5-aminolevulinic acid or an ALA ester, or methyl aminolevulinate (not supplied with the FLUXMEDICARE® system V2

- Product contraindicated for pregnant women.
- Product contraindicated for patients with a history of abnormal photosensitivity.
- Patients allergic to the cream containing 5-aminolevulinic acid or an ALA ester or methyl aminolevulinate or the excipients in the cream (not supplied), including ground-nut oil, peanuts, or soya.
- Porphyria
- Patient who is immunosuppressed for therapeutic reasons, an idiopathic or specific illness

In addition, it is advisable to refer to the package leaflet of the medicines for the specific contra-indications of these specialities.

## 1.5. Undesirable effects

---

The following potential side effects may occur during treatment, they are signs of treatment effectiveness, if in doubt consult your doctor or dermatologist (scabbing, erythema, skin peeling, skin haemorrhage, headache, skin oedema, paresthesia, pruritus, burning sensation of the skin, hot sensation, swelling)

Any serious incident occurring in connection with the FLUXMEDICARE medical device must be reported to MDB TEXINOV and to the competent authority of the Member State in which the user and/or the patient is established.

## 1.6. Precautions for use

---

The FLUXMEDICARE® V2 system is intended for use in combination with a pharmaceutical cream containing 5-aminolevulinic acid or ALA ester, or methyl aminolevulinate (not supplied with the FLUXMEDICARE® V2 system). Please only use creams because other pharmaceutical forms of the active ingredients have not been tested.

The untreated healthy skin beside the lesion does not need to be protected during exposure.

The FLUXMEDICARE® V2 system must be used with a transparent occlusive dressing (not supplied with the FLUXMEDICARE® V2 system). This should carry the CE marking for contraindications and adverse reactions please refer to the instructions for use of this medical device.

The treatment should only be carried out by an authorised medical professional trained in photodynamic therapy.

## 2. Warnings, symbols

---

### 2.1. Warnings and safety remarks

---

This manual contains several warnings describing the banning of or the obligation for certain actions required on the part of the user in order to operate the system in complete safety.

In this manual the warnings are notified by a pictogram  (exclamation mark within a yellow triangle). **Non-adherence to these directions for use may lead to serious injury, to the patient and/or the user!**

In this manual the information is notified by a pictogram  (the letter "i" contained in a blue circle).

**The following warnings should be read carefully!**

#### CAUTION:

- Do not use the device before reading the user manual carefully and understanding it.
- Do not power the device outside the limits specified in this manual or on the label, on the back of the device.
- 
- Risk of eye injuries. Check that safety instructions and principles are followed, and there is no possibility of direct light transmission outside the application area.
- Do not look directly at the light emitted.
- Any opening of equipment requiring a tool is only authorised for MDB Texinov personnel.
- Modifying the FLUXMEDICARE® V2 system is PROHIBITED.
- Do not use for indications other than those specified in this document.
- Risk of eye injuries. Do not leave the light fibres illuminated when these are not being applied to the patient.
- Do not connect accessories other than those specified.
- Do not obstruct the air vents, or place the appliance on a soft or sloping surface.
- Adhere to a free space or a distance from the wall of at least 30 cm on all ventilation faces.
- Do not use an electrical extension lead.
- Only use the count key supplied.

- Do not connect another appliance or count connector (USB type) on the front face of the light source.
- Do not connect this key to a computer or any other electronic appliance.
- Do not place the source on an unstable surface.
- In order to prevent the risk of electric shock, the FLUXMEDICARE® V2 device must be connected to an electric socket with earthing protection.
- The source must not come into contact with a liquid.
- The FLUXMEDICARE® V2 system must not be exposed to high pressures or pointed objects.
- The current socket plug is the disconnection device and should remain accessible.
- Do not use beyond the duration for use.
- In the treatment room, avoid storing and using flammable materials such as oxygen and some anaesthetics.

## 2.2. Safety linked to the light device

Any fault deactivates the light emission.

If the connector disconnects, the light stops, the source emits a sound signal and the treatment is stopped.


## 2.3. Product labelling

Symbol	Meaning
	Manufacturing date
	European Declaration of Conformity
	Batch number
	Reference
	Medical Device
	Manufacturer
	Type B part applied
 	Consult the instruction for use
	Consult the precautions for use
	Laser class 1

Symbol	Meaning
	Temperature limit
	Compliant with directives relating to Waste Electrical and Electronic Equipment (WEEE) and to the limitation of certain hazardous substances in electrical and electronic equipment (RoHS). This symbol affixed to the product or its packaging indicates that the product cannot be discarded in ordinary household waste.
	Expiry date
	ON button (turning on)
	OFF button (turning off)
	Single use: do not re-use the net


## 2.4. Laser warning label

The following label is located on the front part of the light source. It indicates the laser class.

Class1 laser emission

## 2.5. Transport and storage markings

Symbol	Meaning
	The arrows point upwards, store and transport the parcel vertically
	Acceptable temperature limits (-25°C to 50°C)

Symbol	Meaning
	Fragile, avoid knocking and jolting
	Store in a dry place
	Acceptable range of ambient pressure
	Acceptable relative humidity range

### 3. Technical characteristics

#### 3.1. Description

The power emitted does not need to be calibrated by the user because it is controlled by the device itself.

#### 3.2. General information


Parameters	Unit	Values	Comments
<b>Wavelength</b>	nm	638	-
<b>Irradiance</b>	W	2,5	Laser diode output
<b>Irradiance</b>	mW/cm <sup>2</sup>	1.3	On the luminous bands: > 1 mW/cm <sup>2</sup> and <7mW/cm <sup>2</sup>
<b>Treatment duration</b>	min	150	± 1 continuously

Description	Classification
<b>Laser type</b>	Diode laser
<b>Laser class</b>	1
<b>Protection against electric shocks</b>	Class I
<b>Protection against harmful penetration of water or particles</b>	IP 20
<b>Part applied: luminous textile band</b>	Type B

#### 3.3. Electrical and optical performances

Parameter	Range	Unit
<b>Input voltage</b>	220-240	VAC
<b>Input frequency</b>	50-60	Hz
<b>Maximum input current</b>	1	A

### 3.4. Physical characteristics of the light source

### 3.5. Ambient conditions

The system is designed to be used in a standard hospital or clinic environment, with the following ambient conditions:

Parameter	Unit	Min	Max	Comment
Temperature for use	°C	+15	+30	-
Humidity for use	% RH	10	85	Environment for use without condensation required
Atmospheric pressure	hPa	500	1060	-

### 3.6. Transport and storage conditions

Parameter	Unit	Min	Max
Transport and storage temperature	°C	-25	+50
Transport and storage humidity	% RH	10	85
Atmospheric pressure	hPa	500	1060

It is necessary to keep the packaging and protective foam in case the electro-medical device has to be returned. After transport or storage the user must ensure the device is installed in its normal conditions of use (see 3.5. Ambient conditions).

The system must not be stored near a heat source or in direct sunlight.  
The luminous bands must be stored flat; the optical fibres must not be folded or crushed.

Avoid exposure to very dusty environments. Avoid storage and use of flammable materials such as oxygen and some anaesthetics in the treatment room.

The intermediate unit may be held on the main unit using the dedicated housing.

### 3.7. Inspection of the system

---

- i** After unpacking the FLUXMEDICARE® V2 device, inspect it and check it has not suffered any damage during transport.

Before use, make sure the device has not been damaged, that the various parts are intact, e.g.

- The various parts are not broken
- The sheath connecting the main unit of the Medical Device and the intermediate unit should not be damaged in any way (particularly cuts, even partial, pinching). The devices for connecting the sheath to the unit must not be loose.
- The seams on the cap are in good condition
- The textile optical fibres are not damaged

If damage is observed, contact MDB TEXINOV® before use.

### 3.8. Preparation of the source before treatment

---

#### Environment

- i** The equipment must be used by health professionals practising in doctors' surgeries, hospitals and clinics. Do not use in conditions where the ambient temperature exceeds 30°C or where there is condensation.

#### Installation

**i** The FLUXMEDICARE® V2 source from MDB TEXINOV® must be installed on a stable flat surface. Particular attention should be paid to checking that the ventilation outlets are not obstructed (see §2.1). Sufficient space between the front and rear panels and the adjacent surface must be maintained, in order to prevent any obstruction of the air flows that could cause internal overheating of the system.

- i** The protective intermediate unit cap must be removed before the textile is connected.

### 3.9. Connection of the optical fibres

---

The connection between the MDB TEXINOV® source and the optical fibres must be made before starting the system. The connector must be connected to the intermediate unit.

## 4. Protocol for use

---

### 4.1. Preparation of the patient's skin

---

Before applying the photosensitising cream the skin should be disinfected and the surface to be treated cleaned.

Using a spatula apply a layer of photosensitising cream (about 1 mm thick) on a radius of 5 to 10 mm of healthy skin around the lesions.

Cover the area to be treated with a transparent dressing.

All this is left as it is for 30 minutes, avoiding exposure to light.


### 4.2. Installation of the FLUXMEDICARE® V2 TEXLIGHT1-V2 system

---

1. Place the luminous bands (sheath), previously protected with the FMPOX protective cover (cf. installation manual), on the patient's head. The luminous bands are put in place with the luminous face above the lesions to be treated and on the transparent dressing.
2. Position the stretch net (NETOX) on the navy blue side of the sheath so that the luminous bands are fitted on the areas to be treated and to prevent the light from escaping.
3. Put the cap on the net. The cap ties must be positioned behind the head. Tie the cap ties behind the head to adjust the hold and prevent any leak of light during the treatment.
4. The system's various layers must remain positioned as they are for a period of 30 minutes before the start of the light emission. After 30 minutes the adjustment of the cap is checked before starting the light emission.

### 4.3. Putting the intermediate unit and harness in place

1. Place the light source on a flat, stable surface, e.g. a table beside the person to be treated.

2. Put the chest harness in place.  
The attaching harness is completely adjustable. The adjustment is made using the 5 adjustment buckles.

- Open the attaching harness using the clip by pressing it.
- Put the harness on like a waistcoat.
- Close the harness again with the clip.
- Tighten or loosen the buckles on the chest and waist to adjust it.
- Remove the intermediate unit from its support. Remove the protective connector cap and place the flat surface of the intermediate unit facing the patient's chest at the height of the harness.

• Insert the harness buckle at a slant in the intermediate unit slot.

• Slide the harness Velcro into the intermediate unit slot. Then insert the Velcro into the buckle and attach it to the adherent part.


- Check the intermediate unit is correctly attached.
3. Connect the count key supplied with the TEXTLIGHT1-V2 luminous bands in its slot on the top of the installed intermediate unit. This must remain connected throughout the treatment.
  4. Insert the metal end (connector) of the TEXTLIGHT1-V2 luminous bands into the slot in the intermediate unit (check that the protective cap has been removed beforehand).
  5. Reposition the textile part on top of the connect to protect it from light emission during the treatment. The patient is ready to be treated.

#### 4.4. Use of the light source


Diagram 2: Front face of the light source

The front face of the light source (see Diagram 1) is made up of:

- A display screen (user interface)
- A "START/PAUSE" button
- A "STOP" button, a connection for the reinforced cable which connects the source to the intermediate unit
- A blue indicator light when the textile is illuminated (ON/OFF)
- A red indicator light which comes if there is a stoppage in the case of an irradiance problem or when the optical connector is disconnected.


Diagram 2: Rear face of the light source

The rear face of the light source (see Diagram 2) is made up of:

- A connector for the power supply cable
- A start button for the source
- A legal notices label
- A housing for fuses

Diagram 3: Right face of the light source

The right face of the light source (see Diagram 3) is made up of:

- A optical output with reinforced sheath
- A storage support for the intermediate unit
- A legal notices label

#### 4.4.1. Powering up the source

1. Connect the light source to the mains electricity.
2. Turn on the source by pressing the bistable button ON (I)/OFF (O) located at the rear of the light source.  
The display screen comes on.

#### 4.4.2. Display screen

---

The display screen is on the front panel. The verification of parameters and the starting of treatment take place through this interface via the buttons located on the front of the appliance.

The screen displays all the important parameters related to the treatment and the device functions.

#### 4.4.3. Initialisation

---

- During the display of the initialisation screen, the system carries out a self-diagnosis and configures the source parameters.
- The software version is specified on the screen when the source starts up.

#### 4.4.4. Starting treatment

---

- After initialisation, the warm-up of the light source laser is carried out. Treatment can be started from 15°C. This step may take several minutes.
- If the optical connector for the TEXTLIGHTX-V2 luminous bands is not connected, the screen tells the user to connect it to the intermediate unit and the red error warning light comes on. Once the connection has been made, press the START/PAUSE button to continue.
- If the count key is not connected, the screen tells the user to connect it to the intermediate unit with a sound signal and the red error warning light on. Once the connection has been made, press the START/PAUSE button to continue.
- Once the initialisation has finished, the screen displays the number of treatments remaining on the count key.

- To start the treatment, press the “START/PAUSE” button.

- The treatment starts. The blue indicator light comes on to show that the treatment is in progress. The remaining treatment time is shown on the screen (in h = hours: min = minutes: s = seconds).
- The luminous band light up.

#### 4.4.5. Pausing

---

- When the laser is emitting, the blue indicator light is on, the screen displays the remaining treatment time.

- The treatment can be stopped for a short while by pressing the START/PAUSE button. A beeps sounds when the button is pressed and the blue indicator light goes off.

- The screen displays “PAUSE” with a time limit of 15 min. The remaining treatment time is shown underneath.

- To continue the treatment, press the “START/PAUSE” button again.
- If the treatment is not restarted before the end of the pause time (15 min), it is automatically cancelled. The screen shows the end of treatment and a sound signal is heard until the appliance is turned off.

If the patient has to move away from the source (pause), the intermediate unit and textile (combined) group must be place on a stable support.

#### 4.4.6. Stopping before the end of treatment

---

When the laser is emitting the screen displays the remaining treatment time and the blue indicator light is on.

- The treatment can be stopped by pressing the “STOP” button. A beeps sounds when the button is pressed.

- The light emission stops and the next screen is displayed.

- If you want to continue the treatment, press the “START/PAUSE” button.
- If you want to stop the treatment, press the “STOP” button again.
- The screen shows the end of treatment.  
A sound signal is heard until the source is turned off by the user with the ON (I)/OFF (O) button located at the back of the light source.

#### 4.5. End of treatment

- If the treatment proceeds until its end, the screen shows the end of treatment.

- The light emission from the luminous textile stops and the blue indicator light goes off. A sound signal is heard until the source is turned over by the user.
- Use the “ON (I)/OFF (O)” button at the back of the appliance to turn the source off.
- Disconnect the textile and the count key from the intermediate unit.
- Remove the cap, the stretch net, the luminous bands and then the transparent dressing. The transparent dressing and the stretch net must be discarded in the appropriate waste bin.
- Clean the area of treated skin with saline solution.
- Remove the intermediate unit from the harness.
- Put the protective cap back on the connector of the intermediate unit.
- Store the intermediate unit in the support provided on the right face of the light source.
- Remove the harness from the patient.

#### 4.6. Installation of other FLUXMEDICARE® V2 accessories

1. Place the luminous textile (sheath) on the head or face of the patient, with the side with the fibre optic band towards the inside.



The transparent dressing is placed on the lesions to be treated.  
The bands are put in place with the illuminated face against the transparent dressing.

- Position the stretch net on the navy blue side of the sheath so that the luminous bands are fitted on the areas to be treated and to prevent the light from being visible to the exterior.

Reference	Photographs
<p><b>TEXTLIGHT2-V2 or TEXTLIGHT3-V2 Unilateral sheath</b></p>	
<p><b>TEXTLIGHT6-V2 Bilateral temple sheath</b></p>	

## 5. Error messages

Error	Description	Action required
E1	 <p>S/N <small>XXXXXX</small> 0/50 medical@texinov.fr</p>	Contact MDB TEXINOV to obtain new luminous bands.
E2		Connect the luminous band ferrules to the intermediate unit of the source and press the "START/PAUSE" button. While the laser is not reconnected, the red warning light is on. In the event of accidental disconnection during treatment, a sound signal is heard in addition to the illuminated red warning light.
E3		Stop the source and wait until it cools. Check that the air openings are not obstructed and that the ambient temperature does not exceed the nominal values. If the problem persists contact Technical Support. When starting next the appliance will display an error. A diagnosis will be performed.
E4		Stop the source and restart it. If the problem persists contact Technical Support. When starting next the appliance will display an error. A diagnosis will be performed.
E5		The count key is not connected. Connect the count key and press "START/PAUSE". While the key is not reconnected a sound signal is heard and the red warning light is on.

Error	Description	Action required
E6		The count key is not the one that was connected at the start of treatment. Reconnect the count key used initially and press "START/PAUSE". While the initial key is not reconnected a sound signal is heard and the red warning light is on.
		If the screen does not turn on or does not display anything do not use the source. Contact MDB Texinov

If the light source overheats, the rearmable thermal circuit breaker will cut the power supply for the Medical Device automatically. The treatment in progress must be restarted. The Medical Device cannot restart whilst the system is not at nominal operating temperature (see 3 .5). The user must ensure that the ventilation is not obstructed. No action is necessary.

## 6. What to do in case of emergency

In case of emergency, the light source must be cut (disconnected) using the mains plug. In this case, the treatment must not be continued to its end and it is essential that the user contacts MDB TEXINOV to make a diagnosis to find out whether the appliance can be started again.

Any serious incident in connection with the device should be notified to the manufacturer and to the competent authority of the Member State in which the user and / or patient is established.

## 7. Replacement of fuses

### CAUTION:

- If it is necessary to replace the fuses, replace them with fuses with equivalent characteristics:

**F 2AH, 250V**

- Before any fuse replacement, remove the power supply cord from the mains socket.
- If the fuses open again, contact technical support.

- Location of the fuse housing

	<ul style="list-style-type: none"> <li>• Opening of the fuse housing</li> </ul>
	<ul style="list-style-type: none"> <li>• Fuse change</li> </ul>
	<ul style="list-style-type: none"> <li>• Replace with fuses with identical characteristics: F 2AH, 250V</li> </ul>

## 8. Cleaning

---

### 8.1. General

---

#### CAUTION:

- Any cleaning must be carried out when the source is disconnected from the mains electricity.
- The luminous bands are disconnected from the source.
- Never immerse a photodynamic therapy component in a liquid.
- The device must not be sterilised by autoclave or gas.
- Do not bleach.
- Do not dry in a machine.
- Do not iron.
- Do not dry clean.
- Only wash with soapy water.

### 8.2. Source

---

- After turning off the source wait 5 minutes before cleaning.
- Disconnect the power supply lead for the source before cleaning (risk of electric shock).
- Never spray product directly onto the surface of the source.
- Wet a soft, clean cloth with water. Clean the outside of the source and the power supply lead with the cloth. Make sure that no liquid enters the source.
- Dampen a second soft, clean cloth with water. Wipe the source.
- Before use make sure that the source is completely dry. Check that the electric lead is not damaged.

### 8.3. Cap

---

The cap can be washed in the machine on the delicate programme at 40°C without fabric softener.

### 8.4. The stretch net

---

After use discard the single use net.

### 8.5. Luminous bands


---

In case of dirty light strips please follow the instructions below:

- Clean with soapy water
- let the material air dry. DO NOT DRY WITH ARTIFICIAL HEAT.
- Wipe the light strips with a dry cloth.

Before each use, check the condition of the light strips. Particular attention should be paid to:

- The state of cleanliness: If, despite the cleaning described above, the light strips are stained or dirty, they must be discarded.
- The condition of the fibres on the luminous part. Light strips must be discarded if they have more than 10 light points (see photo below)

## 9. Control of the source

---

The source must be checked before each use for obvious visible damage, damaged cables, cracks, worn label.

The source is guaranteed for 2 years against material or manufacturing defects, from the delivery date. During this period, the defective equipment will be replaced free of charge, unless the fault is due to abuse or incorrect use.

**A control of the source must be performed annually** by MDB Texinov® personnel or distributor approved by MDB Texinov.

The service life of the light source is 5000 hours.


MDB TEXINOV® will only take responsibility for the safety, reliability and perfect operation of the appliance if:

- All the adjustments as well as modifications or repairs are carried out by MDB TEXINOV®.
- The device is used in compliance with the instructions for use supplied.

### **Important indications**

Only MDB TEXINOV® personnel are authorised to repair FLUXMEDICARE® V2 products.  
The units must not be opened.

If the equipment must be sent back, it must be packed in its original packaging (see diagram below).

## 10. Recycling

---

Electrical appliances must not be discarded with household waste.

In France, it must be taken to an appropriate collection centre for recycling of electrical and electronic equipment.

To ensure correct treatment of product waste, refer to the usual recycling network, in compliance with local requirements for electrical equipment.

## 11. Electromagnetic compatibility

---

- The use of this appliance beside other appliances should be avoided. This could cause a malfunction. If this use is necessary this appliance, and the others, should be monitored to ensure they are functioning normally.

- The use of accessories, transducers and cables other than those specified or supplied, may cause an increase in electromagnetic emissions or a decrease in the immunity of this machine and cause inappropriate operation.

- RF portable communications devices (including peripherals such as antenna cables and external antennae) should not be used closer than 30 cm (12 inches) from any part of the APPLIANCE, including cables specified by the manufacturer. Otherwise the performance of these devices could be affected.

Electromagnetic compatibility (EMC) and electrostatic discharge (ESD) tests have been performed in accordance with the following standards:

Phenomenon	Standard base used	Result
<b>Radiated emission</b>	CISPR 11 +A1 : 2010	Meets the requirements of Group 1 Class B
<b>Emission by conduction</b>	CISPR 11 Group 1 Class B	Meets the requirements of Group 1 Class B
<b>Harmonic current emission 50 Hz</b>	EN/IEC 61000-3-2+A1+A2 : 2009	Meets the requirements of Class B
<b>Voltage fluctuation</b>	EN/IEC 61000-3-3 : 2013	Meets the requirements
<b>Electrostatic discharges</b>	IEC 61000-4-2 : 2008	Meets the requirements of Class B
<b>Radiated electrostatic fields</b>	EN 61000-4-2 : 2008	Meets the requirements of Class B
<b>Fast transient electrical disturbance/burst</b>	EN/IEC 61000-4-4-3+ A1+A2 : 2010	Meets the requirements of Class B
<b>Shock wave</b>	IEC 61000-4-5 : 2005	Meets the requirements of Class B
<b>RF conducted disturbances</b>	EN 61000-4-6 : 2013	Meets the requirements of Class B
<b>Magnetic field 50 Hz</b>	IEC 61000-4-8 : 2009	Meets the requirements of Class B
<b>Supply voltage interruption and dip</b>	EN/IEC 61000-4-11 :2004	Meets the requirements of Class B

## 12. Regulation

The laser source:

- is a Class IIa appliance compliant with the regulations on medical products
- is a laser Class 1 laser device (in accordance with the standard, NF EN 60825-1),

has CE marking in accordance with regulation 2017/745.

The FLUXMEDICARE® V2 device is a medical device CE marked in June 2020 by the BSI organisation (CE 2797) in accordance with the regulation 2017/745 and is in risk class IIa.

The accessories are CE marked Class I medical devices in accordance with Regulation (EU) 2017/745.

The following standards are applicable.

Reference	Title
<b>NF EN ISO 13485</b> <b>version: September 2016</b>	Medical Devices - Quality Management Systems - Requirements for Regulatory Purposes
<b>Regulation (EU) 2017/745</b>	New European regulation of 05 April 2017
<b>NF EN 60825-1</b> <b>version: October 2014</b>	Safety of Laser Products - Part 1: Equipment classification, requirements and user's guide
<b>ISO 15223-1</b> <b>version: ISO 15223-1</b> <b>July 2021</b>	Symbols to be used with medical device labels, labelling and information to be supplied
<b>NF EN 60601-1/A1</b> <b>version: March 2014</b>	Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance

Reference	Title
<b>CEI 60601-1-2:2014</b> version: February 2014	Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
<b>NF EN 60601-1-6</b> version: July 2010	Medical Electrical Equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability
<b>NF EN 62304</b> version: October 2006	Medical device software - Software life-cycle processes
<b>NF EN 62304/A1</b> version: April 2018	Medical device software - Software life-cycle processes Amendment A1
<b>NF EN 62366-1 + AMD1</b> version July 2020 version: December 2015	Medical devices - Part 1: Application of usability engineering to medical devices
<b>NF EN ISO 14971</b> version: December 2019	Medical devices - Application of risk management to medical devices

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